

Hand-Activated Monopolar Insulated Forceps



Fine tips for accuracy and speed of tissue dissection



Hand-activated switch eliminates the need for foot pedals



Attached cord can be used with most ESUs

Our fine-point, accurate, electrosurgical tissue forceps have an enhanced design with even more sophistication.

This is the new standard instrument for both tissue dissection and hemostasis – saving time and money.

It's now available in two lengths for complete body and facial plastic surgery procedures.

B 65244	10"	25.5 cm
B 65245	7"	17.5 cm



5238 Royal Woods Parkway
Tucker, GA 30084
877.252.2517
770.414.4880
770.414.4879 Fax
www.blackandblacksurgical.com

To order or for additional information, contact Black & Black customer service at 877.252.2517 or 770.414.4880 or visit our on-line web store at www.blackandblacksurgical.com.